## PAS - Template - Record Keeping

## Steriliser Monitoring

Use either a Class 4, 5 or 6 chemical indicator, or record observations for <u>every cycle</u>. Temperature observations are to be recorded every 15 seconds during key stages of Type B or S steam sterilisation cycles, and at least every 30 seconds during key stages of a displacement steam or dry heat sterilisation cycle. Refer to the content of AS4815:2006 Appendix G.

Date	Cycle number or start time	Cycle type (B, S etc)	Product/s (empty or list items)	Monitoring by observation (record every 15-30 seconds)				Monitoring by Indicator			
				Time	Temp °C	Pressure kPa	Result	Indicator type (Class 4/5/6)	Indicator Batch/Lot Number	Result	Operator
										1	
										1	